

Bath | circulation water bath


circulation water bath, 6 / 11 / 22l, up to 100°C ±0.1°C, 5 l/min, flat lid included

IQ | OQ
available
page 12 for details

WCB-22 with flat lid (included)

Stainless steel flat lid (included)

WCB 100°C, 5 l/min, 6/11/22 l

Ideal for:

- biotechnology, clinical, environmental, medical, petroleum & food industry, pharmaceutical and industrial applications
- connection with evaporators and viscometers is applicable

Features:

- user's self-compensation function to control difference between real temperature and indicated temperature (±10.0°C)
- stainless steel bath & lid for superior durability & high thermal efficiency / consistency
- storage function for the set values of temperature and timer
- powerful circulation pump assures temperature uniformity: internal and external circulation
- storage function for the set values of temperature and timer
- backlight LCD display
- CE certified and unique serial number for tracing

Safety mechanism:

- locking mode
- over-temperature and over-current protection
- sensor error detection
- low water level warning and automatic stop function



Controller:

- digital fuzzy controller with **Jog-Shuttle-switch** (turn + push)

Model	WCB-6	WCB-11	WCB-22
Capacity	6 l	11 l	22 l
Temperature range & accuracy	room temperature +5°C - 100°C, ±0.1°C		
Temperature resolution, uniformity & probe	0.1°C-display, 0.1°C-control, ±0.2°C, PT100		
Heating power	800 W	1 kW	2 kW
Pump	max. flow rate 5 l/min (internal/external)		
Circulation nozzle	inner Ø 7 x outer Ø 11 mm (inner Ø 9/outer Ø 12 mm silicon tube connection for inner circulation)		
Timer & alarm	99hr 59 min (delay & continuous run), error status & timer end		
Drain valve	-		mounted (Ø9.5 - 11 mm)
Effective space	120 x 155 x 150 mm	170 x 240 x 150 mm	372 x 295 x 150 mm
Internal dimensions (W x D x H)	300 x 155 x 150 mm	302 x 240 x 150 mm	500 x 295 x 150 mm
External dimensions (W x D x H)	330 x 185 x 385 mm	335 x 275 x 385 mm	540 x 335 x 410 mm
Packing size & gross weight	350 x 210 x 340 mm, 11 kg	360 x 300 x 420 mm, 13 kg	580 x 380 x 470 mm, 18 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz		
Order number 230V	DH.WCB00106	DH.WCB00111	DH.WCB00122
Order number 120V	DH.WCB00206	DH.WCB00211	DH.WCB00222

Bath | circulation water bath

circulation water bath, 8 / 12 / 22 / 30 l, -25°C - +150°C ±0.1°C, 25 l/min, flat lid included

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- connection with evaporator

Features:

- stainless steel bath, polyurethane insulation of 40 mm thickness and body made of powder-coated steel
- innovated refrigeration system ensures temperature performance, includes intelligent control for preventing compressor overload, CFC-free (R-404A)
- innovated powerful circulation pump ensures temperature uniformity: internal and external circulation and 5 steps circulation speed
- stainless steel flat lid & valve ball included, dome lid optional
- RS232C interface: remote monitoring and controlling with PC
- proper liquid level tag in bath
- user's self-compensation function to control difference between real temperature and indicated temperature (±10.0°C)
- backlight LCD display
- CE certified and unique serial number for tracing

Safety mechanism:

- overheat and over-current protection
- sensor error detection and leakage breaker
- low liquid level warning and automatic stop function
- locking mode



Controller:

- digital fuzzy controller with **Jog-Dial-switch** (turn + push) & capacitive **touch button**

Model	WCR-8	WCR-12	WCR-22	WCR-30
Capacity	8 l	12 l	22 l	30 l
Temperature range & accuracy	-25°C - +150°C, ±0.1°C			
Temperature resolution	0.1°C-display, 0.1°C-control			
Temperature uniformity & probe	±0.2°C at -10°C, PT100			
Heating Power	1.0 kW	1.5 kW	2 kW	2.2 kW
External circulation nozzle	inner Ø10.5 mm, outer Ø12.5 mm, for inner circulating: connect directly the in-/outlet nozzles			
Pump	max. flow rate: 25l/min, max. pressure: 7.5 psi, max. lift: 4 m			
Refrigerator	250 W	250 W	370 W	560 W
Cooling capacity at +20°C	300 W	380 W	680 W	850 W
Cooling capacity at 0°C	200 W	250 W	500 W	620 W
Cooling capacity at -20°C	110 W	130 W	300 W	400 W
Timer & alarm	99hr 59 min (delay & continuous run), error status & timer end			
Effective space	175 x 153 x 150 mm	230 x 181 x 150 mm	290 x 220 x 200 mm	330 x 290 x 200 mm
Internal dimensions (W x D x H, mm)	175 x 320 x 150 mm	230 x 350 x 150 mm	290 x 390 x 200 mm	330 x 460 x 200 mm
External dimensions (W x D x H, mm)	330 x 474 x 752 mm	400 x 523 x 807 mm	483 x 563 x 850 mm	494 x 633 x 906 mm
Packing size & gross weight	560 x 420 x 922 mm, 45 kg	630 x 530 x 972 mm, 54 kg	670 x 615 x 1022 mm, 63 kg	740 x 625 x 1072 mm, 70 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 230V	DH.WCR00408	DH.WCR00412	DH.WCR00422	DH.WCR00430
Order number 120V	DH.WCR00508	DH.WCR00512	DH.WCR00522	DH.WCR00530

Bath | circulation water bath

circulation water bath, 8 / 12 / 22 / 30 l, -35°C - +150°C ±0.1°C, 25 l/min, flat lid included

IQ | OQ
available
page 12 for details

WCL -35°C - +150°C, digital

Ideal for:

- connection with evaporator

Features:

- stainless steel bath, polyurethane insulation of 40 mm thickness and body made of powder-coated steel
- innovated refrigeration system ensures temperature performance, includes intelligent control for preventing compressor overload, CFC-free (R-404A)
- innovated powerful circulation pump ensures temperature uniformity: internal and external circulation and 5 steps circulation speed
- stainless steel flat lid & valve ball included, dome lid optional
- RS232C interface: remote monitoring and controlling with PC
- proper liquid level tag in bath
- user's self-compensation function to control difference between real temperature and indicated temperature (±10.0°C)
- backlight LCD display
- CE certified and unique serial number for tracing

Safety mechanism:

- overheat and over-current protection
- sensor error detection and leakage breaker
- low liquid level warning and automatic stop function
- locking mode

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Contents

1 Introduction	13
2 Safety Instructions	13
2.1 Description of symbols	13
2.2 Instructions.....	13
3 Product description	15
3.1 Key features.....	15
3.2 General survey	15
3.3 Operating surface	16
3.4 Term description	16
4 Packing content	17
5 Positioning the device and commissioning	17
6 Job description	18
6.1 Menu structure	18
6.2 Set turn-off delay.....	18
6.3 Set turn-on delay.....	18
6.4 Set temperature.....	18
6.5 Justify.....	18
6.6 Set pump power	19
6.7 Key lock.....	19
6.8 Storage function.....	19
7 Maintenance and cleaning	20
8 Troubleshooting Guide	20
9 Technical data	21
10 Certificate	22

1 Introduction

Thank you for purchasing our **WCH High Temperature Bath Circulator**. You did a good choice buying a **Wisdom** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a safe place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.

- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.
- Before you switch on the device, make sure that the connections (see 3.2) for external water circulations are connected with each other or with an external device.

- It is important that the device stands up on a horizontal, stabile and firm subsurface to avoid vibrations.
- Make sure that the device can not slip away.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.
- Make sure the device stands upright within the last 24 hours.
- Do not tilt the device.

- Do not use this device in area with potentially explosive atmosphere.

- Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted up to 30°C over max. temperature. Otherwise the heating maybe doesn't work.

- Do not put heat-sensitive items at the device.
- Warn all attendees if you leave the device unattended.
- Pay attention when you lift up the lid, hot condensation water may drop down or hot steam could come out.
- Do not clean the scorching device.


- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.
- Before you switch on the device, make sure that the connections (see 3.2) for external water circulations are connected with each other.

3 Product description

3.1 Key features

- Simplest operating by MAIN button and Jog-Dial knob.
- Fast commissioning.
- Long life by simply and solid construction.
- Highest temperature accuracy by digital fuzzy controller up to 250 °C.
- Temperature and delay time showed on illuminated LCD.
- Storage function for set time and temperature values.
- Turn-off an -on delay from one minute up to over 99 hours in one minute steps adjustable.
- Signal sound after end of turn-off delay.
- Easy to justify.
- Key lock against inadvertent changing of set speed and timer values.
- Over temperature protection by safety switch.
- A lot of laboratory devices are mountable.

3.2 General survey

3.3 Operating surface

3.4 Term description

SV: set value (upper line, in key lock menu no SV)

PV: present value, measured by the device (lower line, in key lock menu also upper line)

Main switch

- On / off switch.

Timer lamp

- Is blinking in menu items for timer modulation.
- Shines continuous in the key lock menu item when timer is activated.

Heating lamp

- Is blinking in menu item for temperature modulation.
- Shines if the heater is activated.

Pump lamp

- Is blinking in menu item for pump modulation.
- Shines continuous in the key lock menu item if pump is on.

Jog-Dial knob

- To change values by turning and stores, confirms, stops and starts by pushing.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot. Thos function is
- Set the maximal temperature roughly 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

MAIN button

- For changing between the different menu items and the key combination.

SUB Button

- To reach the sub menu Items.

Connections for external water circulation

- To connect external devices.   If you directly connect both connections the media is circulating inside the device.

4 Packing content

quantity	denomination
1	WCH High Temperature Bath Circulator
1	power supply cord
1	stainless steel lid
1	silicon tube
1	operating manual


5 Positioning the device and commissioning

1. Remove the package and perhaps soiling carefully.
2. If there are any unexpected transport damages, please inform your forwarder and your customer service.
3. Set the protection switch 30 °C over desired bath temperature otherwise the heating coil does not work.
4. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatized.
5. Connect the connections for external water circulation with each other or connect external water cycle.
6. Make sure the external water cycle is proof.
7. Make sure that you follow all safety instructions.
8. Connect the device with the power outlet.
9. Fill the operating liquid in.
10. Switch the device on and start the process (look at job description).

6 Job description

6.1 Menu structure

6.2 Set turn-off delay

The turn-off delay timer serves as automatically stop function. It turns off the heating and pump automatically.

- First menu item after turning on is “set turn-off delay”. The timer lamp is blinking.
- Turn the Jog-Dial to set the desired turn-off time.
- Start and stop timer by pushing the Jog-Dial knob.
- Timer must be stopped to enable a renewed changing of timer value.
- In the upper display line is the SV (set value) and in the lower line the remaining time to turn-off (PV) shown.
- Stop timer for endurance run.

6.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating and pump. It is possible to work with turn-off timer in addition.

- Have a look to the menu overview (6.1) to open the relevant menu item.
- Turn the Jog-Dial to set the desired turn-on time.
- Start and stop timer by pushing the Jog-Dial knob.
- Timer must be stopped to enable a renewed changing of timer value.
- In the upper display line is the SV (set value) and in the lower line the remaining time to turn-on (PV) shown.
- Stop timer for immediately run.

6.4 Set temperature

- To reach this menu item you have to push the MAIN button in the menu “set turn-off delay”
- The temperature is adjustable in 0,1°C-steps.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the

- Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.
 - Turn the Dial knob clockwise to increase the temperature and counter-clockwise to reduce it.
 - The heating lamp shines, if the heater is activated.

6.5 Justify

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about $\pm 10^{\circ}\text{C}$.

The Offset is calculated as follows:

Temperature measured by you – showed Temperature (PV) = Offset

- To reach this menu you have to push the SUB button in the menu item “set temperature”.
- You can recognize this menu at the PV showed in the SV-line.
- If you turn the Jog-Dial knob.
- Set the calculated value by turning the Jog-Dial. The Offset will be added to the PV immediately.
- Save and confirm the Offset by pushing the MAIN button, you will also reach the previous menu.

6.6 Set pump power

- To reach this menu item you have to push the MAIN button in the menu “set temperature”.
- To set pump power turn the Jog Dial knob. The performance is adjustable in steps of 1-5. To confirm you have to press the Jog Dial knob.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.

6.7 Key lock

- To avoid accidental changing of the set values, push the MAIN button in the menu item “set temperature”.
- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay in the PV-line. The function of the Jog-Dial knob is deactivated.
- Press the MAIN button once again to leave the key lock. Than you are in the menu item “set turn-off delay”.

6.8 Storage function

Storage Function is very useful if you want to work with the same parameters often. Then you didn't need to set them in every usage.

- For storing values you only need to hold the Jog-Dial knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jug-Dial till “0” and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired value and use the Jog-Dial knob for storing.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	<ul style="list-style-type: none"> ✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is “ON”. ✓ Fuse is blown.
Unit Not Heating	<ul style="list-style-type: none"> ✓ Check if the MAIN Power Switch is “ON”. ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 6.4) ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).
Heating, but unit not reaching temperature	<ul style="list-style-type: none"> ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV). ✓ Problem in the sensor or main controller.
ERROR 1	<ul style="list-style-type: none"> ✓ Temperature Sensor Error.
ERROR 4	<ul style="list-style-type: none"> ✓ Error in circulation.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by [witeg](#) authorized people. Contraventions exclude the warranty.

9 Technical data

		WCH-8	WCH-12	WCH-22	WCH-30
Dimension (B×T×H)	Bath (mm) Effective space	180x130x155	250x130x155	340x250x155	340x250x200
	Bath (mm)	180x320x155	250x320x155	340x440x155	340x440x200
	External (mm)	303x443x569	373x443x569	463x563x569	463x563x614
Capacity		8 Lit.	12 Lit.	22 Lit.	30 Lit.
Heating power		1 kW	1,5 kW	2,2 kW	2,5 kW
Temp.	Range	Ambient + 5 °C to + 250 °C			
	Accuracy	± 0,1 °C			
	Uniformity	± 0,2 °C			
	Sensor	PT100			
Pump		Max Flow Rate : 16 Lit/min, Max Pressure : 5psi Pump Max. Height: 2,8 m			
Controller		Digital fuzzy controller with Jog-Dial knob (Turn + Push)			
Display		Digital LCD with Back-Light Function			
Timer		99 h 59 min			
Safety device		Over temp. and over current protection, low filling level detector			
Material	Internal	Stainless steel			
	Capsule	Powder coated steel			
	Lid	Stainless steel			
Others		<ul style="list-style-type: none"> - Storage function (temperature and timer), - Key lock (Jog-Dial input disabled) - Alarm (sensor error, low filling level and timer-end) 			
Packing Size (mm)		410x470x630	410x470x630	500x580x690	500x580x690
Gross Weight		19 kg	21 kg	23 kg	37 kg
Power supply		AC 230 V, 1N~, 50/60 Hz			

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	WCH High Temperature Circulation Water Bath
Typ	
Serialno.	
Date	


witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de

Bath | circulation water bath

circulation water bath, 6 / 11 / 22 / 30 l, up to 250°C ±0.1°C, 16 l/min, flat lid included

WCH 250°C, 16 l/min, 6/12/22/30 l

Ideal for:

- applications at high temperature
- biotechnology, clinical, environmental, medical, petroleum & food industry, pharmaceutical and industrial applications
- connection with evaporators and viscometers is applicable

Features:



- corrosion resistant stainless steel bath with 1.2 mm thickness
- ceramic fiber insulator of 50 mm thickness and powder-coated steel body
- user's self-compensation function to control difference between real temperature and indicated temperature (±10.0°C)
- digital timer function (99hr 59min): delay time and operation time can be controlled
- drain valve mounted and stainless steel flat lid included
- storage function for the set values of temperature and timer
- backlight LCD display
- CE certified and unique serial number for tracing

Safety mechanism:


- overheat and over-current protection
- sensor error detection
- locking mode

Controller:

- digital fuzzy controller with **Jog-Shuttle-switch** (turn + push)

Stainless steel flat lid (included)


Connection for external circulation Ø12.5mm

Model	WCH-8	WCH-12	WCH-22	WCH-30
Capacity	8 l	12 l	22 l	30 l
Temperature range & accuracy	room temperature +5°C - 250°C, ±0.1°C			
Temperature resolution	0.1°C-display, 0.5°C-control			
Temperature uniformity & probe	±0.2°C, PT100			
Heating Power	1.0 kW	1.5 kW	2.2 kW	2.5 kW
External circulation nozzle	inner Ø10.5 mm, outer Ø12.5 mm			
Pump	max. flow rate: 16l/min, max. pressure: 5 psi, max. lift: 2.8m			
Timer & alarm	99hr 59 min (delay & continuous run), error status & timer end			
Effective space	180 x 130 x 155 mm	250 x 130 x 155 mm	340 x 250 x 155 mm	340 x 250 x 200 mm
Internal dimensions (W x D x H)	180 x 320 x 155 mm	250 x 320 x 155 mm	340 x 440 x 155 mm	340 x 440 x 200 mm
External dimensions (W x D x H)	303 x 443 x 569 mm	373 x 443 x 569 mm	463 x 563 x 569 mm	436 x 563 x 614 mm
Packing size & gross weight	470 x 400 x 500 mm, 19 kg	470 x 400 x 500 mm, 21 kg	600 x 490 x 500 mm, 23 kg	600 x 490 x 500 mm, 37 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 230V	DH.WCH00408	DH.WCH00412	DH.WCH00422	DH.WCH00430
Order number 120V	DH.WCH00508	DH.WCH00512	DH.WCH00522	DH.WCH00530

Bath | circulation water bath

circulation water bath, 6 / 11 / 22 / 30 l, up to 250°C ±0.1°C, 16 l/min, flat lid included

WCH 250°C, 16 l/min, 6/12/22/30 l

Ideal for:

- applications at high temperature
- biotechnology, clinical, environmental, medical, petroleum & food industry, pharmaceutical and industrial applications
- connection with evaporators and viscometers is applicable

Features:


- corrosion resistant stainless steel bath with 1.2 mm thickness
- ceramic fiber insulator of 50 mm thickness and powder-coated steel body
- user's self-compensation function to control difference between real temperature and indicated temperature (±10.0°C)
- digital timer function (99hr 59min): delay time and operation time can be controlled
- drain valve mounted and stainless steel flat lid included
- storage function for the set values of temperature and timer
- backlight LCD display
- CE certified and unique serial number for tracing

Safety mechanism:

- overheat and over-current protection
- sensor error detection
- locking mode

Controller:

- digital fuzzy controller with **Jog-Shuttle-switch** (turn + push)

Stainless steel flat lid (included)

Connection for external circulation Ø12.5mm

Model	WCH-8	WCH-12	WCH-22	WCH-30
Capacity	8 l	12 l	22 l	30 l
Temperature range & accuracy	room temperature +5°C - 250°C, ±0.1°C			
Temperature resolution	0.1°C-display, 0.5°C-control			
Temperature uniformity & probe	±0.2°C, PT100			
Heating Power	1.0 kW	1.5 kW	2.2 kW	2.5 kW
External circulation nozzle	inner Ø10.5 mm, outer Ø12.5 mm			
Pump	max. flow rate: 16l/min, max. pressure: 5 psi, max. lift: 2.8m			
Timer & alarm	99hr 59 min (delay & continuous run), error status & timer end			
Effective space	180 x 130 x 155 mm	250 x 130 x 155 mm	340 x 250 x 155 mm	340 x 250 x 200 mm
Internal dimensions (W x D x H)	180 x 320 x 155 mm	250 x 320 x 155 mm	340 x 440 x 155 mm	340 x 440 x 200 mm
External dimensions (W x D x H)	303 x 443 x 569 mm	373 x 443 x 569 mm	463 x 563 x 569 mm	436 x 563 x 614 mm
Packing size & gross weight	470 x 400 x 500 mm, 19 kg	470 x 400 x 500 mm, 21 kg	600 x 490 x 500 mm, 23 kg	600 x 490 x 500 mm, 37 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 230V	DH.WCH00408	DH.WCH00412	DH.WCH00422	DH.WCH00430
Order number 120V	DH.WCH00508	DH.WCH00512	DH.WCH00522	DH.WCH00530