



Heating mantle for flasks “Fibroman-D”

FEATURES

Soft mantle made with a double woven mineral wool. A mains switch and safety earth are also fitted. The mantle can be supported by a retort ring or tripod stand.

MODELS

| Part No. | Flask capacity ml | Ø / Height (exterior) cm | | Power W | Weight Kg |
|----------|-------------------|--------------------------|-----|---------|-----------|
| 3004192 | 250 | 15 | 6.5 | 130 | 0.8 |
| 3004195 | 500 | 17 | 7 | 270 | 0.9 |
| 3041910 | 1000 | 20 | 8.5 | 410 | 1 |
| 3041920 | 2000 | 24 | 11 | 530 | 1.5 |

See accessories (page 230) and controller (pages 293 and 294).

ACCESSORIES

| Retort Part No. | Tripod Ring Part No. | Mantle Part No. |
|-----------------|----------------------|-----------------|
| 7041921 | 7041922 | 3004192 |
| 7041951 | 7041952 | 3004195 |
| 7419101 | 7419102 | 3041910 |
| 7419201 | 7419202 | 3041920 |



Heating mantle for flasks and funnels “Fibroman-O”

WITH TRIPOD SUPPORT STAND.

APPLICATIONS

Suitable for round bottom flasks, decanting funnels and reaction flasks with a bottom drain.

FEATURES

- Woven mineral fibre mantle.
- Epoxy coated exterior.
- On/Off switch.
- Safety earth.

MODELS

| Part No. | Flask capacity ml | Ø / Height (Exterior dims.) cm | | Ø cm Bottom hole | Tripod Height cm | Power W | Weight Kg |
|----------|-------------------|--------------------------------|------|------------------|------------------|---------|-----------|
| 3049110 | 1000 | 24 | 10.5 | 8 | 34 | 410 | 1.5 |
| 3049120 | 2000 | 28.5 | 12 | 8 | 36 | 530 | 2.5 |
| 3049150 | 5000 | 37.5 | 17.5 | 9 | 38 | 840 | 3.5 |
| 3491100 | 10000 | 43.5 | 20 | 9 | 43 | 1400 | 5 |
| 3491200 | 20000 | 48 | 25 | 12 | 36 | 2300 | 13 |

Due to its capacity the part no. 3491200, 20000 ml capacity is manufactured in a cylindrical case. See accessories (page 230) and controller (pages 293 and 294).



Heating mantle for flasks

ADAPTABLE FOR THE “AGIMAN” STIRRER UNIT.

MODELS

| Part No. | Flask capacity ml | Ø / Height (Exterior dims.) cm | | Power W | Weight Kg |
|----------|-------------------|--------------------------------|-----|---------|-----------|
| 3000447 | 250 | 15 | 6.5 | 130 | 0.8 |
| 3000448 | 500 | 17 | 7 | 270 | 0.9 |
| 3000449 | 1000 | 20 | 8.5 | 410 | 1 |
| 3000450 | 2000 | 24 | 11 | 530 | 1.5 |

ACCESSORIES

Mantle magnetic stirrer “Agiman”. With interchangeable mantle adapters. Electronic control of stir speed (up to 1600 r.p.m.). Part No. **7000446** See page 26. With connection at the back for an electronic contact thermometer “Sensoterm” or digital electronic regulator “Electemp-TFT”. See controller (pages 293 and 294).



Stirring mantle with adapter device.



Stirring mantle, complete.