



# HOTPLATES



## Round Hotplates

HOTPLATE TEMPERATURE UP TO 400 °C. WITH 6 PLACE HEATING POWER CONTROLLER.

### FEATURES

Cast steel plate with embedded heating elements throughout the plate surface.

Fast heating rate at maximum power setting.

Enamelled outer casing that is resistant to high temperatures and corrosive agents.

### CONTROL PANEL

Controller 6 position power settings.



### MODELS

Part No.	Plate Ø cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
1000718	12	8 20 20	800	2.1
1000442	15	8 22 22	1600	3
1000443	19	8 26 26	2100	4.5