

0.1 mg - 220 g



WBA-220 0.1mg - 220g analytical

Ideal for:

- analytical balancing

Features:

- various counting & weighing modes (g, ct, etc.)
- with Ø80 mm/Ø90 mm weighing plate and glass draft shield
- weighing response time adjustable
- electromagnetic base
- numeric linearity calibration
- built-in RS232C port for printer
- blue digital LCD display with backlight function
- 200g calibration weight included (only with Ext-CAL model)
- stabilization time (typical at fast) ≤ 8 seconds
- CE certified and unique serial number for tracing



200g calibration weight included WBA-220

Model	WBA-220 Ext-CAL	WBA-220A Int-CAL
Weighting range, weighting plate	0.1 mg - 220 g, Ø80 mm	0.1 mg - 220 g, Ø90 mm
Read out, reproducibility & linearity	0.1 mg	
Permissible temperature	5.0 - 35°C	
Weight modes	g / ct	g / ct / % / GN / dwt / etc.
Counting mode	counting & weighing mode available	
Calibration	external calibration, inkl. 200g Kalibergewicht	internal calibration, AUTO
Dimension & weight	190 x 240 x 265 mm, 6.3 kg	223 x 345 x 331 mm, 5.7 kg
Packaging size	290 x 420 x 490 mm, 8.3 kg	390 x 510 x 520 mm, 10.1 kg
Power supply	DC 9V with AC 230V adapter, 50/60 Hz	DC 15V with AC 230V adapter, 50/60 Hz
Order number	DH.WBA0220	DH.WBA0220A

NEW



Thermal printer, for balance WBA, RS232 connection for paper rolls (1 paper roll included)

Article	Order number
WBA-100 printer	DH.WBA100P
5x spare paper roll for WBA-100	DH.WBA0084



Thermal printer for balance WBA, RS232 connection for label paper rolls (1 label paper roll included)

Article	Order number
WBA-420 printer	DH.WBA420P
5x spare label paper roll for WBA-420	DH.WBA0771




Sample plates, 10 pieces

Diameter	Order number
Ø90mm	DH.WBA0005

Balance

high-precision lab balance, 0.001g - 320/640g, percent weighing mode

0.001 g - 320 g
- 640g



WBA-620

WBA -320 /-620 /-620A high-precision

Ideal for:

- high precision balancing

Features:

- percent (%) weighing mode available
- with Ø90 mm/Ø118 mm/Ø128 mm weighing plate and glass draft shield (included)
- weighing response time adjustable
- electromagnetic base
- numeric linearity calibration
- built-in RS232C & USB port for printer
- super-size backlight LCD display, can display present temperature
- time and date setting for recording data
- stabilization time (typical at fast) ≤ 8 seconds
- battery operation possible (6V, 1.2A)
- CE certified and unique serial number for tracing

Safety mechanism:

- four level anti-shock function adjustable
- default and auto lock present working mode



200g (WBA-320) / 500g (WBA-620)
calibration weight included

NEW

Model	WBA-320 Ext-CAL	WBA-620 Ext-CAL	WBA-620A Int-CAL
Weighting range, weighing plate	0.001 - 320 g, Ø 90 mm	0.001 - 620 g, Ø 128 mm	0.001 - 620 g, Ø 118 mm
Read out, reproducibility & linearity	0.001 g		
Permissible temperature	5.0 - 35°C		
Weight modes	% mode / g / ct / GN / dwt / etc.		
Counting mode	available counting & weighing mode		
Calibration	external calibration, 200g weight included	external calibration, 500g weight included	internal calibration, AUTO
Dimension & weight	223 x 345 x 331 mm, 3.6 kg	223 x 345 x 331 mm, 5.7 kg	
Packaging size	290 x 385 x 365 mm, 5.1 kg	390 x 510 x 520 mm, 10.1 kg	
Power supply	DC 12V with AC 230V adapter, 50/60 Hz	DC 12V with AC 230V adapter, 50/60 Hz	
Order number	DH.WBA0320	DH.WBA0620	DH.WBA0620A



Thermal printer, for balance WBA, RS232 connection
for paper rolls (1 paper roll included)

Article	Order number
WBA-100 printer	DH.WBA100P
5x spare paper roll for WBA-100	DH.WBA0084



Thermal printer for balance WBA, RS232 connection
for label paper rolls (1 label paper roll included)

Article	Order number
WBA-420 printer	DH.WBA420P
5x spare label paper roll for WBA-420	DH.WBA0771




Sample plates, 10 pieces

Diameter	Order number
Ø90mm	DH.WBA0005

Balance

standard lab balance, 0.01g - 3200g/6200g, big weighing plate, various counting & weighing modes



WBA 0.01g - 3200 g/6200 g standard



WBA-3200

Ideal for:

- standard balancing purposes

Features:

- various counting & weighing modes (g, ct, GN, dwt, etc.)
- with 168 x 168 mm weighing plate and optional glass draft shield
- weighing response time adjustable
- electromagnetic base
- external calibration and numeric linearity calibration
- built-in RS232C & USB port for printer
- super-size backlight LCD display, can display present temperature
- stabilization time (typical at fast) ≤ 8 seconds
- 2 kg calibration weight included (only WBA-3200)
- CE-MDD certified and unique serial number for tracing

Safety mechanism:

- four level anti-shock adjustable function
- default and auto lock present working mode



2 kg (WBA-3200) calibration weight included

Model	WBA-3200	WBA-6200
Weighting range, weighing plate	0.01 - 3200 g, 168 x 168 mm	0.01 - 6200 g, 168 x 168 mm
Read out, reproducibility & linearity	0.01 g	
Permissible temperature	5.0 - 35°C	
Weight modes	% mode/ g / ct / GN / dwt /etc.	
Counting mode	available counting & weighing mode	
Calibration	external calibration, 2 kg calibration weight included	external calibration, <u>no</u> weight included
Dimension & weight	195 x 295 x 90 mm, 2.6 kg	
Packaging size	243 x 377 x 163 x mm, 3.2 kg	
Power supply	DC 12V with AC 230V adapter, 50/60 Hz	
Order number	DH.WBA3200	DH.WBA6200



Thermal printer, for balance WBA, RS232 connection for paper rolls (1 paper roll included)

Article	Order number
WBA-100 printer	DH.WBA100P
5x spare paper roll for WBA-100	DH.WBA0084



Thermal printer for balance WBA, RS232 connection for label paper rolls (1 label paper roll included)

Article	Order number
WBA-420 printer	DH.WBA420P
5x spare label paper roll for WBA-420	DH.WBA0771



Glass draft shield for WBA-3200/6200 three slide doors (left, right, top)

Order number
175 x 195 x h230 mm
DH.WBA0003

Balances

Accessories: thermal printer, sample plates, glass draft shield

Thermal printer, for balance WBA, RS232 connection
for paper rolls (1 paper roll included)

Article	Order number
WBA-100 printer	DH.WBA100P
5x spare paper roll for WBA-100	DH.WBA0084



Thermal printer for balance WBA, RS232 connection
for label paper rolls (1 label paper roll included)

Article	Order number
WBA-420 printer	DH.WBA420P
5x spare label paper roll for WBA-420	DH.WBA0771



Sample plates, 10 pieces

Diameter	Order number
Ø90mm	DH.WBA0005



Glass draft shield for WBA-3200/6200
three slide doors (left, right, top)

	Order number
175 x 195 x h230 mm	DH.WBA0003

