

Galaxy 330 with 4 display panels

Ideal for:

- total bacteria counting
- total coliform counting


Features:

- 4 displays showing present count and three previous counts
- possibility of average value
- storage for up to 100 counts
- subtraction function
- suitable for petri dishes from Ø90 to Ø150 mm and Ø60 mm
- pressure sensor allows to operate with any marking pen
- CE certified



Galaxy 330

Model	Galaxy 330
Display	3-digit (0-999) LED display, 4 panels
Magnification	2 ~ 3 X
Petri dish diameter	Ø90-150 mm adjustable, Ø60 mm
Memory	up to 100 counts
Average	Yes
Zero reset	Yes
Count black	Yes
Illumination	LED (adjustable)
Pressure sensitivity	Adjustable
Acoustical counting signal	Adjustable
Output	USB
Dimension (LxWxH)	35.1 x 30.2 x 37 cm
Net weight	5.6 kg
Order number EU plug	00829 175330-02
Order number US plug	00829 175330-01



Model	Software
Purpose	Add, delete, search, print data as well as save as csv file for usage in a spreadsheet
Scope of supply	Software including USB stick and connection cable
Order number	00829 175330-12