

# Cooling Water Circulator

## HS-3005N (CE), HS-3005, HS-3005W

As a source of cooling water to the Rotary Evaporator and as a heat-emitting equipment it supplies sufficient quantity of cooling water as highly efficient chemical pumps are installed and does not contaminate other instruments as it does not get corroded.



MODELS	HS-3005N (CE)	HS-3005	HS-3005W
Capacity	4 Liter	15 Liter	24 Liter
Bath Dimension (ØD × H)	190 × 200 mm	280 × 270 mm	350 × 250 mm
Refrigerator	1/3 HP	3/4 HP	2 HP
Dimension (W × D × H)	250 x 510 x 510 mm	330 × 560 × 670 mm	630 × 680 × 900 mm
Gross Weight	28 Kg	40 Kg	110 Kg
Power Consumption	334 W	750 W	2,000 W
Circulator Pump	Magnet 20 W Max. flow 4.5 L/min		Magnet 40 W max. flow 15 L/min
Refrigerant	134 A		R 404
Temp. Range	Down to ~ -20 °C		
Controller	FND Controller		
Power Supply	AC 230 V, 50/60 Hz		