

Rotator

5 - 60 rpm, 0 - 90° angle, multiple discs available

RT-10 5 - 60 rpm, 0 - 90° angle



RT-10 with disk WRT120

Ideal for:

- general mixing, diffusing, dialyzing, hybridization and aerating culture
- keeping biological samples in suspension

Features:

- 5 - 60 rpm
- 0-90° angle adjustable
- digital LCD display
- smooth start and change of stirring speed while operation
- timer for 0-99 hours
- 5 different disks available
- durable construction for long life time
- CE certified and unique serial number for tracing

Safety mechanism:

- timer 99h 59min with timer-end signal
- operation LEDs
- overload protection
- locking mode

Controller:

- programmable **Jog-Shuttle** controller with touch-button

Model	RT-10
Design	turning motion
Speed range	5-60 rpm
Angle range	0-90°
Programs	6 programs, 10 steps
Timer & alarm	99hr 59min, error status timer-end
Display	digital LCD with backlight function
Controller	digital feedback controller (Jog-Shuttle)
Temperature resistance	lower than 65°C (dry air)
Safety circuit	overload protection
Dimension & weight	200x250xh260 mm, 6 kg
Packaging size	370x240xh430 mm, 8.5 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz
Order number 230V	DH.WRT00010
Order number 120V	DH.WRT00020

Please note: no disks included - please order separately!

Scope of supply:

- Rotator RT-10
- 4x knobs, 4x posts WRT300



4x knobs, 4x posts for installing disks (included in scope of supply)


Article	Order number
WRT300	DH.WRT00300

Accessories:



Rotator disc for 8x 50ml conical tubes

Article	Order number
WRT100	DH.WRT00100



Rotator disc for 16x 15ml conical tubes

Article	Order number
WRT110	DH.WRT00110



Rotator disc for 32x 1.5ml conical/centrifuge tubes

Article	Order number
WRT120	DH.WRT00120



Rotator disc for 60x 1.5ml conical/centrifuge tubes

Article	Order number
WRT130	DH.WRT00130



Rotator disc for 4x 250ml flasks

Article	Order number
WRT140	DH.WRT00140