

Vortex mixer

up to 3,300 rpm, orbital shaking motion, multiple heads available

VM-10
with platform head



VM-10 orbital motion

Ideal for:

- mixing solutions in test tubes or small flasks with circular, oscillating movement

Features:

- 0 - 3,300 rpm
- orbital motion
- inclusive platform head Ø76 mm and pop-off cup head
- 3-position power switch for continuous run or „touch-on“
- stepless adjustable speed by turning the dial-knob
- tough metal housing and reliable suction cups ensure stable motion, low noise and safe working
- various shapes, sizes and materials of the head allow mixing of almost all common tubes and containers
- CE certified and unique serial number for tracing

Controller:

- **analog phase** controller
- 3-position switch for start, stop and „touch-on“
- dial-knob for stepless speed setting

Model	VM-10
Design	vortex mixer
Speed range	0-3,300 rpm
Shaking motion	orbital
Controller	analog phase controller
Material	aluminium die-casting with powder-coating
Head	included pop-off cup head and Ø76 mm platform
Dimension & weight	160x128x145 mm, 3 kg
Packaging size	260x200x220 mm, 4 kg
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz
Order number 230V	DH.WVM00010
Order number 120V	DH.WVM00020

Scope of supply:

- Vortex mixer VM-10
- platform head PM110 Ø76mm
- pop-off cup head PM210
- removing tool for heads



Platform head PM110



Removing tool



Pop-off cup head PM210

Accessories:



Microtube platform head for 0.5 - 2.0 ml

Article	Order number
PM310	DH.WVM00310