

Compact Shaker KS 15 A



Description

The reliable Orbital Shaker KS 15 A for general shaking tasks for daily use. The loading capacity is 15 kg on a shaking platform of 400 x 300 mm.

Basic Equipment

Basic device without shaking plate, **without rack system**

Loading Capacity

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	15 pieces
Erlenmeyer flasks 500 ml	12 pieces

Erlenmeyer flasks 1000 ml 5 pieces
Separating funnels 250 ml 1 pieces

You will find further loading possibilities [here](#)

Details

Motion: orbital
Shaking platform: 400 x 300 mm
Max. load: 15 kg
Shaking speed: 30 - 420 rpm, steplessly adjustabl
Stroke: 17 mm
Timer: 0 - 120 min / continuous
Order No.: 6171 000
Accessory: Rack System Combifix KS
Universal tray KS
Incubator Hood TH 15

Technical specifications

Electrical supply: 230 V or 115 V, 50 / 60 H please indicate in case of order
Enclosure protection: IP 21
Heat emission: approx. 5 W
Ambient temperature: 5°C to 50°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 510 x 490 x 150 mm
Weight: 18 kg