

## Compact Shaker KS 15 B



### Description

The reliable Horizontal Shaker KS 15 B for general shaking tasks for daily use. The loading capacity is 15 kg on a shaking platform of 400 x 300 mm.

### Basic Equipment

Basic device without shaking plate, **without rack system**

### Loading Capacity

|                          |           |
|--------------------------|-----------|
| Erlenmeyer flasks 100 ml | 20 pieces |
| Erlenmeyer flasks 250 ml | 15 pieces |

Erlenmeyer flasks 500 ml 12 pieces  
Erlenmeyer flasks 1000 ml 5 pieces  
Separating funnels 250 ml 1 pieces

You will find further loading possibilities [here](#)

## Details

**Motion:** reciprocating (horizontal)  
**Shaking platform:** 400 x 300 mm  
**Max. load:** 15 kg  
**Shaking speed:** 30 - 420 rpm, steplessly adjustable  
**Stroke:** 17 mm  
**Timer:** 0 - 120 min / continuous  
**Order No.:** 6173 000  
**Accessory:** Rack System Combifix KS  
Universal tray KS  
Incubator Hood TH 15

## Technical specifications

**Electrical supply:** 230 V or 115 V, 50 / 60 H please indicate in case of order  
**Enclosure protection:** IP 21  
**Heat emission:** approx. 5 W  
**Ambient temperature:** 5°C to 50°C  
**Relative humidity:** ~ 85 %  
**Dimensions (w x d x h):** 510 x 490 x 150 mm  
**Weight:** 18 kg