

Dual-Action Shaker KL 2



Description

Depending on the task the shaking motion of the Dual Action Shaker KL 2 can be switched from orbital to reciprocating.

Basic Equipment

Basic device incl. rack system with 3 clamping strips and rubber mat

Loading Capacity

Erlenmeyer flasks 250 ml 9 pieces

Separating funnels 250 ml 1 pieces

You will find further loading possibilities [here](#)

Details

Motion: can be switched from orbital to reciprocating (horizontal)

Shaking platform: 300 x 300 mm

Max. load: 2 kg

Shaking speed: 30 - 420 rpm, steplessly adjustable

Stroke: 8 mm

Timer: 0 - 120 min / continuous

Order No.: 6115 000

Accessory: Clamping strip h 0050 118

Clamping strip v 0050 477

Cramp strip 0050 206

Spring strip 0050 207

Test tube rack with 44 drillings à 14 mm Ø

Test tube rack with 44 drillings à 16 mm Ø

Test tube rack with 44 drillings à 18 mm Ø

Test tube rack with 14 drillings à 30 mm Ø

Technical specifications

Electrical supply: 230 V or 115 V, 50 / 60 Hz please indicate in case of order

Enclosure protection: IP 21

Heat emission: approx. 4 W

Ambient temperature: 5°C bis 50°C

Relative humidity: ~ 85 %

Dimensions (w x d x h): 355 x 455 x 195 mm

Weight: 13,5 kg