

Three-Dimensional Shaker TL 2



Description

The Laboratory Shaker TL 2 works with slow, three-dimensional motion for gentle, uniform wetting. Well suited for petri dishes, culture bottles etc.

Basic Equipment

Basic device incl. basic tray and rubber mat

Details

Motion: three-dimensional orbital motion

Shaking platform: 220 x 220 mm

Max. load: 2 kg
Shaking speed: 1 - 22 rpm, steplessly adjustable
Angle of tilt: 8° fixed
runtime: 0 - 9 h 59 min / continuous
Order No.: 6139 000
Accessory: Additional Tray TL 2

Technical specifications

Electrical supply: 100 - 240 V, 50 / 60 Hz
Enclosure protection: IP 21
Heat emission: max. 0,5 - 2 W
Ambient temperature:: 5°C to 45°C
Relative humidity: ~ 85 %
Dimensions (w x d x h): 250 x 415 x 118 mm
Weight: 9,1 kg