

## Tilt Shaker WS 10



### Description

The Tilt Shaker WS 10 works with slow, linear tilting motion and a fix angle of tilt of 5° with a shaking speed of 1 - 23 rpm. The shaking platform of 340 x 340 mm can be loaded with a weight of 10 kg.

### Basic Equipment

Basic device incl. basic tray and rubber mat

## Details

<b>Motion:</b>	linear tilting motion
<b>Shaking platform:</b>	340 x 340 mm
<b>Max. load:</b>	10 kg
<b>Shaking speed:</b>	1 - 23 rpm, steplessly adjustable
<b>Angle of tilt:</b>	5° fixed
<b>Timer:</b>	0 - 120 min / continuous
<b>Order No.:</b>	6114 000
<b>Accessory:</b>	Multi-storey rack system WS / TL

## Technical specifications

<b>Electrical supply:</b>	24 V DC, 50 / 60 Hz 115 V version on request
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 3 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	355 x 450 x 205 mm
<b>Weight:</b>	18,5 kg